

Stand-up aid with accessory adaptor to suit Sentida SC

Code: 20043

The stand-up aid has been designed with the safety of residents in mind. The loop style stand-up aid fixes to the side of the bed frame in multiple positions providing stability when getting up from the bed. It supports residents to move independently, limiting physical strain and preventing the risk of falls. The stand-up aid includes an accessory adaptor to fit the Sentida SC bed.

FEATURES

- Includes accessory adaptor to suit Sentida SC bed
- Provides safety and stability when exiting the bed
- Limits physical strain for residents
- Suits multiple positions



SPECIFICATIONS

Code(s)	20043
Overall dimensions	W 300 H 512mm
Weight	1.37kg
Safe working load	75kg (top), 50kg (side)

Ready to order?

Our customer support team are available 5 days a week to help answer all your questions.

Phone: 0800 656 527

Phone: +64 7 578 7228