

Mira 90 procedure light

Mobile



Amico introduces the latest advancement in minor procedure and diagnostic lighting technology with the Mira 90 procedure light. Precisely positioned LED's allow for superior shadow control, depth of field and light intensity. A CRI of 95 and colour temperature of 4500K produce excellent colour accuracy for minor procedures and examinations.


Being a true cross-platform design, the Mira 90 procedure light can be used in your ED, Trauma, ICU, LDR and procedure room environments.

Model	Code
Amico Mira 90 procedurelight - Mobile	0320-1990-05
Amico Mira 90 procedure light - Wall mounted	0320-1993-05
Amico Mira 90 procedure light - Ceiling mounted single	0320-1996-05
Amico Mira 90 procedure light - Ceiling mounted dual	0320-1998-05

Wall mounted



Ceiling mounted (single)



Ceiling mounted (dual)



Features

Light intensity of 90,000 lux at 1 metre

Twist handle for focusing light beam from 16.51cm to 30.5cm

Control pad on light head for easy access to light functions

Slim design to reduce clutter in small workspaces

Lightweight long reach arms allow for easy positioning

Removable/sterilisable handle

Specifications

Voltage	120 – 277 volts
Shadow Dilution	Single mask: 0.7%,
	Cavity: 100%
	Single Mask with Cavity: 0.3%
	Double Mask with Cavity: 55.3%
Light head diameter	45.72cm
Light field diameter	16.51cm - 30.5cm
Weight of lighthead	5kg
Brightness control	50-100%
Depth of illumination	240.03cm
Light Intensity (Central luminance)	90,000 lux, at 1m
Electrical	100-240 VAC, 50/60 Hz, 0.42 amps
Color Rendering Index (Ra)	95
Standard colour temperature	4500K
Number of LEDs	36
Life of LEDs	50,000 hours
Mounting options	Ceiling, wall, floor stand
Certifications	cETLus conform to UL60601-1; UL60601-2-41; UL certified IP54 per IEC 60598
Warranty	Five years



Made in Canada

Get in touch & connect   

Email sales@cubro.co.nz • Call us 0800 101 255 • Visit cubro.co.nz

Cubro[®]

© Copyright 2023 Cubro Ltd.