

# DeVilbiss 10 litre compact oxygen concentrator

Code: 813

The drive DeVilbiss 10L oxygen concentrator offers optimal oxygen therapy for those requiring varying flow rates. Designed with enhanced quality, durability and simplicity in mind, the patented DeVilbiss Turn-Down Technology minimises wear on internal components and reduces power consumption by at least 15%, extending the life of the concentrator. It utilises the latest advancements in oxygen therapy and is compatible to be used with cylinder transfill systems.

## FEATURES

- Small, lightweight and compact design, weighing only 19kg
- Capable of delivering 87% to 96% of O<sub>2</sub> purity from 2 to 10 LPM
- Oxygen outlet incorporating a fire protection adapter for added safety
- Visual and audible alarms to indicate power failure, high gas temperature, high pressure, low/high flow, low oxygen and service required
- An auxiliary oxygen port is included which can be used to fill cylinder transfill systems
- Durable and reliable, tested under high temperatures and high humidity
- Bright LED panel for patient controls, with a protected cannula fitting and recessed humidifier nook to prevent damage
- Lockable flowmeter that can accommodate up to 15m of oxygen tubing plus 7ft of cannula



## Ready to order?

Our customer support team are available 5 days a week to help answer all your questions.

Phone: 0800 656 527

Phone: +64 7 578 7228

Email: [sales@cubro.co.nz](mailto:sales@cubro.co.nz)

## DeVilbiss 10 litre compact oxygen concentrator

### SPECIFICATIONS

Code(s)	813
Outlet pressure	20.0 ± 1.0 psi (138 kPa +/- 7 kPa)
Oxygen percentage	2-10 LPM=87%-96%
Power consumption	120 vac, 60HZ: 639 watts average
Weight	19kg
Dimensions	H 622 W 342 D 304mm
Warranty	3Years
Alarms	Audible and visual high pressure, low/high flow, low oxygen, high gas temperature, power failure, Oxygen Sensing Device (OSD)